

PROFESSIONAL LIGHT TECHNOLOGY • CLINICAL VERSATILITY • BUSINESS VALUE

RADIANT GENESIS

5-IN-1 DPL PROFESSIONAL LIGHT PLATFORM

ONE DEVICE. FIVE IN-DEMAND TREATMENT PATHWAYS.

A professional, Health Canada-licensed Class III platform designed to give trained practitioners greater control over light delivery, thermal management and client comfort.



VASCULAR



FIRMING



PIGMENT



HAIR



ACNE



PRECISION IS THE DIFFERENCE.

Built for clinics and medical-aesthetic practices that want more than a basic, single-purpose IPL system.



01 • PROFESSIONAL POSITIONING

Not a basic IPL. A precision-minded professional platform.

Radiant Genesis is designed to elevate light-based treatments from a generic flash of energy into a controlled professional service. Its dual application handpieces, adjustable sub-pulse architecture, large treatment window and active cooling system give the practitioner meaningful control over how energy is selected, sequenced and managed at the skin surface.



VERIFIED CANADIAN MEDICAL-DEVICE STATUS

DPL versus conventional IPL

IPL describes a broadband, polychromatic light source. Filters and settings tailor that light toward targets such as melanin and haemoglobin. The professional value of Radiant Genesis is not that it is simply 'stronger'; its hardware and pulse controls are designed to make delivery more purposeful, repeatable and comfortable.

CLINICAL QUESTION	CONVENTIONAL / BASIC IPL	RADIANT GENESIS DPL PLATFORM
 How is light applied?	Broad-spectrum flashes shaped by filters and general settings.	Two purpose-specific handpieces: HR 640–1200 nm and SR 530–1200 nm.
 How is energy timed?	Pulse options vary by model and may be limited.	Pulse-string output with adjustable 1–20 ms sub-pulses and 1–4 s intervals.
 How is heat managed?	Cooling capability varies widely between systems.	Air, water and semiconductor cooling; handle temperature 0–15°C.
 How quickly can areas be covered?	Spot size and repetition rate depend on the platform.	A large 10 × 50 mm treatment area supports efficient coverage.
 What is the practical advantage?	Generic settings can limit customization.	More ways to match handpiece, pulse and fluence to the treatment goal.



TECHNICAL ACCURACY BUILDS TRUST



02 • TREATMENT VERSATILITY

Five clinical pathways. One intelligent platform.

Radiant Genesis gives a clinic the ability to build a broad light-based treatment menu around five of today's most requested aesthetic concerns. Each pathway begins with professional assessment, appropriate settings and a treatment plan tailored to the client.



01 VISIBLE VEINS & REDNESS

Targeted light settings may help reduce the appearance of superficial vascular concerns and persistent redness.



02 SKIN FIRMING & REJUVENATION

Light-based protocols can support collagen-focused services for texture, firmness and a refreshed appearance.



03 PIGMENTATION & SUN DAMAGE

Purposeful wavelength selection can be aligned with visible sun spots and uneven tone when clinically appropriate.



04 LONG-TERM HAIR REDUCTION

The HR handpiece supports professional hair-reduction protocols across appropriate body areas and client profiles.



05 ACNE-FOCUSED SKIN TREATMENTS

Acne-focused light services can be incorporated into a broader professional skin-management plan.



ONE CONSULTATION CAN REVEAL MORE THAN ONE OPPORTUNITY



03 • CLINICAL CONTROL

Designed for diverse skin presentations— not one-size-fits-all settings.

Radiant Genesis is promoted for use across a wide range of skin tones, including darker skin tones. The professional advantage is its capacity for customization: trained operators can adjust pulse width, energy and delivery while using active cooling to support a more controlled treatment experience.

No light-based platform makes every treatment suitable for every person. Fitzpatrick type, melanin content, recent sun exposure, medication, target depth, anatomical site and treatment objective must guide candidacy and settings. This is where a high-level device—and high-level training—matter most.

Why parameter control matters

- Pulse width helps shape how rapidly energy is delivered to the target.
- Fluence determines the energy density applied per square centimetre.
- Pulse sequencing can build therapeutic heat while allowing controlled intervals.
- Cooling supports comfort and surface thermal management.
- Handpiece selection aligns the wavelength band with the intended application.



CONTROL + COOLING + TRAINING

Verified Radiant Genesis architecture

SYSTEM FEATURE	PUBLISHED SPECIFICATION
Professional interface	15-inch true-colour touchscreen
Application handpieces	HR 640–1200 nm • SR 530–1200 nm • 2 handles
Treatment window	10 × 50 mm
Pulse architecture	1–20 ms sub-pulse width • pulse-string output • 1–4 s interval
Energy density	3–25 J/cm ² in 1 J/cm ² steps
Thermal management	Air + water + semiconductor cooling • handle cooling 0–15°C

 **THE RIGHT PROMISE: CONTROLLED, CUSTOMIZED, PROFESSIONALLY ASSESSED.**



04 • PRACTICE VALUE

The value is not one treatment. It is the menu the device unlocks.

A single-purpose device limits how often a treatment room can earn. Radiant Genesis supports five distinct service pillars, multiple body areas, series-based protocols and natural cross-service conversations—all from one professional platform.



**MORE
REASONS
TO BOOK**



**MORE
WAYS
TO RETAIN**



**MORE
VALUE
PER ROOM**

Illustrative service-menu value

Radiant's treatment menu demonstrates the breadth of billable services clinics may build around the platform. The examples below are suggested Canadian menu prices—not revenue guarantees—and should be adapted to the local market, treatment time, practitioner credentials and scope.

SERVICE PILLAR	EXAMPLE SERVICE	SUGGESTED PRICE
Skin rejuvenation	Face / face + neck / face + neck + chest	\$199 / \$249 / \$349
Hair reduction	Upper lip through full legs	\$49–\$299
Vascular appearance	Nose, cheeks and small-to-large areas	\$99–\$399
Pigmentation	Single spot, hands or full face	\$99–\$199
Acne-focused	Face, chest or back	\$149–\$249

The real return on a premium platform

- **Protocol range:** A broader menu reduces reliance on one seasonal service.
- **Client experience:** Active cooling can support comfort and confidence.
- **Workflow:** A large spot size can improve coverage efficiency on appropriate areas.
- **Professional credibility:** Verified Canadian licensing and published specifications support a stronger sales conversation.
- **Growth potential:** Series, maintenance plans and complementary skincare can deepen lifetime client value.



05 • READY-TO-USE POSITIONING

Radiant Genesis: precision light, elevated.

A PROFESSIONAL PLATFORM BUILT TO DO MORE

Professional-use and claims guidance

This overview is sales and educational material, not a treatment protocol. Use only within applicable scope of practice, current device labelling and Health Canada-licensed indications. Follow manufacturer training, contraindication screening, informed consent, eye protection, patch testing and documentation requirements. Results and number of sessions vary. Avoid unqualified claims such as 'safe for everyone,' 'works on every skin type,' 'permanent removal,' 'cures acne' or 'guaranteed results.'



“

RECOMMENDED CLAIM LANGUAGE

“Radiant Genesis is a professional DPL/IPL platform with adjustable pulse controls, two wavelength-specific handpieces, active cooling and a large treatment window. Treatment suitability, settings, comfort and outcomes depend on professional assessment, current device labelling and operator training.”

Selected sources & claim notes

1. Health Canada Medical Devices Active Licence Listing. Radiant Genesis, licence 113227, Class III, identifier RAD100; listing accessed August 6, 2026.
2. Schoenewolf, Barysch & Dummer. 'Intense pulsed light.' Current Problems in Dermatology 42 (2011): 166–172.
3. Radiant Genesis DPL IPL Brochure. Published wavelength bands, pulse settings, treatment window and cooling system; accessed August 6, 2026.
4. Zhao et al. DPL terminology in peer-reviewed literature; cited for terminology only, not as a Radiant Genesis device study.
5. Radiant Genesis DPL Treatment Menu & Pricing Guide. Illustrative Canadian service-menu prices, not income projections.



BOOK A PROFESSIONAL DEMONSTRATION

1-866-330-9490 • INFO@RADIANTLASER.CA • [RADIANTLASER.CA](https://www.RADIANTLASER.CA)

